

## Slit lamp – adaptations

DS-STD is adaptable to common slit lamps in the market:

Manufacturer	Model(Magnification of oculars)	Camera position	Increase of slit lamp length
66 Vision	YZ5S (12.5x)	Downwards	48.6mm
	YZ5F1 (12.5x)		
	YZ5G (12.5x)		
	YZ5X1 (12.5x) Zeiss type		
Argus	SL-5000 (12.5x)	downwards	48.6mm
Bobes	2070 (15x)	downwards	48.6mm
BON	SL-75 ELITE (12.5x)	downwards	48.6mm
	SL-75 (12.5x)		
	SL-85 (12.5x)		
Burton	860(15x)	downwards	48.6mm
CSO	SL-980 (12.5x)	downwards	48.6mm
	SL-990		
HAI	SL-5000 (12.5x)	downwards	48.6mm
Huvitz	HS-5000 (12.5x)/HIS-5000 (12.5x)	downwards	48.6mm
Inami	Classicus L-0179 (15x)	downwards	48.6mm
	Classicus L-0179 Z (15x)		
	Classicus L-0219 (12.5x)		
	L-0187 (15x)		
	L-0189 (12.5x)		
	L-0229 (12.5x)		
	L-0240 (15x)		
SL-0229-E (12.5x)			
Keeler	SL-40H (12.5)	downwards	48.6mm
Marco	Ultra G5 (16x)/Ultra G5 Zoom (16x)	downwards	48.6mm
Nidek	SL-1600 (10x/16x)	downwards	48.6mm
	SL-1800 (12.5x)		
Mediworks	S350S (12.5x)	downwards	48.6mm
	S350C (12.5x)		
	S350 (12.5x)		
	S360 (12.5x)		
	S360S (12.5x)		
Nikon	FS-3 (12.5x)	downwards	48.6mm
	FS-2 (12.5x)		
Reichert	Xcel-255 (12.5x)	downwards	48.6mm
	Xcel-200/250/400/700 (12.5x)		
Righton	RS-1000 (12.5x)	downwards	48.6mm
Rodenstock	RO-2000 (10x)	downwards	48.6mm
	RO-3000 (12.5x)	downwards	48.6mm
	RO-4000 (12.5x)		
	RO-5000 (16x)		
Shin Nippon	SL-203	downwards	48.6mm
Takagi	SM-70 (16x)	downwards	48.6mm
	SM-90 (12.5x)		

Tomey	TSL-3000 (15x)	downwards	48.6mm
Topcon	SL-D2/D4/D4Z/D5 (12.5x)	downwards	48.6mm
	SL-5D (12.5x)/6E (12.5x)	downwards	48.6mm
Woodlyn	HR3 (15x)	downwards	48.6mm
	HR5 (12.5x)		
Zeiss	30 SL (12.5x/10x)	downwards	48.6mm
	SL-130 (12.5x/10x)	downwards	48.6mm
	SL-160 (12.5x/10x)		

Please note, the DS-STD is **NOT** adaptable to the following slit lamps (RED color is HS type)

Manufacturer	Model	Manufacturer	Model
66 Vision	YZ5E (交角式) YZ5J (交角式) YZ5F (交角式) YZ5X (交角式)	Nikon	NS-1
BON	SL-72 SL-2E SL-Zoom	Reichert	Prelude Model 12521 Prelude Model 12522
Burton	800 830 835	Rodenstock	RO-1000
CSO	SL-901/R 980 Zoom / 990 Zoom	Shin Nippon	SL-45 SL-102
Haag-Streit	BM-900 BM-900V BC-900 BD-900 BQ 900 (12.5x)	Takagi	SM-10N / -30N SM-2 (16x) SM-4 (16x)
HAI	SL-5000	Tomey	TSL-1000 TSL-2000
Inami	Classicus L-0169 L-0185 L-0391 (ATTO) L-0392(ATTO) (12.5x) L-0393(ATTO) (15x)	Topcon	SL-1E SL-3D SL-2D/ZE (12.5x) SL-2E (12.5x) SL-2F/D7 (12.5x)
Marco	2b 5 Ultra G2 (16x) Ultra G4 (16x)	Zeiss	10 SL SL 115 20 SL (12.5x / 10x) SL 100/16 (16x) SL 125/16 (16x) SL-105 (12.5x / 10x) SL-120 (12.5x / 10x)
Mediworks	S260S/S260 S280/S280S/S280C	Woodlyn	HR1
Nidek	SL-250 SL-450	Zeiss Jena, produced in GDR	RSL-110 RSL-180
Mentor	BM-900 Version		